


YQC-2L Emergency bed



- It has a separable surface with the main trolley frame which could be removed from trolley and used as a separate stretcher, easy locking/unlocking system to lock the separable stretcher into the trolley
- Swing-down side rails to protect patient when transfer
- Adjustable height by crank
- With IV pole
- $\phi 150\text{mm}$ castors, two with breaks
- It is made of high-strength aluminum alloy
- The trolley is used for carrying patients in the hospitals, emergency centres.

Optional:

Adjustable backrest



YQC-2R Emergency bed

- Swing-down side rails to protect patient when transfer
- Adjustable height by crank
- Adjustable backrest by gas spring
- With IV pole, mattress
- $\phi 150\text{mm}$ centre control castors
- It is made of steel and ABS plastic
- The trolley is used for carrying patients in the hospitals, emergency centres.

