







# YDC-4A/4B

## Scoop Stretcher




### Technical features:

- The stretcher can be divided into left and right halves by the clutches in both ends
- It has extra low profile for sliding the stretcher easily beneath the patient to minimize the risk of further injury.
- Adjustable length according to patient's height, easy to lock and unlock, and foldable
- It is made of high-strength aluminum alloy.
- Can be used together with head immobilize
- Mainly used for ambulance, hospitals, emergency center
- Accessories: 3pcs belts packed in a small bag



# YDC-4D Scoop Stretcher



## *Technical features:*

- the stretcher can be divided into left and right halves by the clutches in both ends
- It has extra low profile for sliding the stretcher easily beneath the patient to minimize the risk of further injury.
- Adjustable length according to patient's height, easy to lock and unlock, and foldable
- It is made of polyethylene ,allow X-ray
- Can bed used together with head immobilize
- Mainly used for ambulance, hospitals, emergency center
- Accessories: 3pcs belts packed in a small bag