


YDC-7A1/7A3/7A4 Spine Board




Technical features:

- It is made of polyethylene material, no contaminator, seamless ,durable, allow X-ray.
- It can be used together with head immobilize.
- It's mainly used for hospitals, sports, ambulance and outdoor carrying patients and wounded person

YDC-7B Spine Board



YDC-7B1



YDC-7B2

Technical features:

- It is made of polyethylene material, no contaminator, seamless ,durable, allow X-ray.
- It can be used together with head immobilize.
- It's mainly used for hospitals, sports, ambulance and outdoor carrying patients and wounded person

Spide Strap

- Secure patients to backboards.
- The size can be adjusted by velcrosuitable for adult and children.
- The width of starp is 5cm.

